



1120 Via Callejon – San Clemente, CA 92673 USA
TEL: +1-949-492-2600 FAX: +1-949-492-1499 e-mail: sbt@southbaytech.com
www.southbaytech.com

Model 38001 TEM Cross Sectioning Kit

38001 Cross Section Kit includes the following:

01-02113-01 Bonding Jig (used for gluing the specimen stack together w/epoxy)	Qty: 1
01-02111-01, 01-02110-01, 01-02112-01 Specimen Holder Pedestal Assembly (used for holding specimen stack during ultrasonic cutting)	Qty: 1
01-02114-01 Brass Reinforcing Tube (used for reinforcing the specimen tube)	Qty: 10
01-02116, 01-02115-01 Teflon Tube Assembly Holder (used for holding the brass tube and specimen tube during gluing together w/epoxy)	Qty: 1
Model 110 Clamp (used for gluing the specimen stack together)	Qty: 1
01-02071-01 Stainless Steel Cutting Tool, 4 x 5 mm (used for cutting out the Si squares for making the specimen stack)	Qty: 1
01-02027-06 Stainless Steel Cutting Tool, 2.3 mm diam. (used for cutting out the specimen tube)	Qty: 1
MWH 135 mounting wax	Qty: Sample