



1120 Via Callejon – San Clemente, CA 92673 USA
TEL: +1-949-492-2600 FAX: +1-949-492-1499 e-mail: sbt@southbaytech.com
www.southbaytech.com

Specifications for SoniCut™ 380 UltraSonic Cutter

Dimensions:	6"W x 19"H x 9"D
Net Weight:	20 lbs.
Transducer Material:	Titanium
Transducer Oscillation:	26 kHz
Tool Material:	Stainless Steel
Frequency control:	SoniLoc™ phase lock loop circuitry
Termination method:	Electronic auto stop or manual termination
Cutting depth control:	Analog dial gauge with 10 micron graduations
Specimen Location Adjustment:	Unlimited with magnetic base plate movement
Front panel controls:	Separate controls for power/illumination and transducer activation
Electrical Input:	100 VAC 50/60Hz 110 VAC 50/60Hz 220 VAC 50/60Hz

Model 36001 Alignment Microscope Specifications (Optional Accessory)

Dimensions:	5"W x 11"H x 8"D
Net Weight:	6 lbs.
Magnification:	40x (other magnifications available)
Field of view:	4mm